

Antigen self-test ("living room test") – Instructions for registering at the Corona Test Portal

Status as of: June 2021

Your pleasant and above all safe stay in Vorarlberg is very important to us. These comprehensive measures are intended to reduce the spread and thus keep the risk of infection to a minimum.

The TVR rule, "Tested, Vaccinated, or Recovered" is valid throughout Austria for admission to accommodations, bars & restaurants, swimming pools, theatres, and many other cultural, sports and recreational facilities.

More information is available at: www.vorarlberg.travel/visiting-safely

During your stay in Austria, you have several options for testing: public testing stations, pharmacies, or simple antigen self-tests. You can easily perform a rapid antigen test ("living room test") yourself and record it in the official Corona Test Portal. You will need the free test kit, a pen, a mobile phone, a tablet or PC with a camera, and an internet connection. You can get the free test kit at the accommodation provider or at the local town hall/tourist office. You will receive the negative test confirmation via the Corona Test Portal. It is valid for 24 hours to provide admission.

Note for foreign guests:

Self-tests registered in the Corona Test Portal are **not valid for border crossing**. For this purpose, an official rapid antigen test or PCR test is required. If needed, reserve a testing appointment well in advance at a public testing station or pharmacy via the Corona Testing Portal at:

<https://covid-portal.lwz-vorarlberg.at>



**Corona Test Hotline:
Tel. 0800 201 360**

1

One-time user registration at the Corona Test Portal

- Register at the state's Corona Test Portal by visiting: <https://covid-portal.lwz-vorarlberg.at>
- Create your personal user account by clicking on the "**Register New User**" button. Next, fill in all necessary data in the registration form, agree to the terms of data usage, and confirm by tapping on "**Register.**"
- You will then receive your user name and a **confirmation link** via email to the email address you provided. Follow the link to activate your login for the Corona Test Portal.
- Remember your login details for additional test registrations at the site.

2

Perform the self-test and register your result

1.

Carry out the self-test according to the instructions for use. Stick the enclosed QR code (valid only once) to the lower area on the test strip so that the test result and QR code are visible. If there is already a QR code on the test strip, you can skip this step.

2.

Open the state's Corona Test Portal by visiting <https://covid-portal.lwz-vorarlberg.at> and then register your self-test. Click on "home self-test" and follow the registration process step by step.

3.

Based on your test result, select "**Result is NEGATIVE**" or "**Result is POSITIVE**"

4.

If the test result is **NEGATIVE**, next activate your camera by clicking on "Start Camera" and capture the test strip within the red frame until it is validated by the system.

5.

Proceed to the next step. Ensure that the test cannot be reused by blacking out or crossing out the test result field. The QR code must remain on the test strip.

6.

Tap on "Start Camera" again and position the self-test within the red frame.

7.

Complete the self-test registration by clicking on "Next Step." **All finished!**

8.

If, contrary to expectations, the test result is **POSITIVE**, make an appointment for an obligatory PCR test directly at the Corona Test Portal, inform your accommodation provider, call the **1450** hotline, and remain in your room.

3

Download the test confirmation

- After registering for the test, you will receive a notification via SMS and email.
- You can **download your test confirmation** directly from the **Corona Test Portal** or from the **email attachment**.
- The test confirmation of your self-test is valid for **24 hours** and serves in printed or digital form as proof for restaurant visits, hotel/accommodation services, as well as access to events or various recreational facilities.